



HEMIGLOBINCYANIDE SOLUTION (HiCN)

FOR USE AS CALIBRATOR WITH THE HEMOGLOBIN DETERMINATION

- Independent Calibrator
- Expiration date > 5 years, after production
- Ready for use
- Liquid



Products	Product no.	Quantity
Hemiglobincyanide Solution (HiCN) Low level	2457	12 x 5 ml
Hemiglobincyanide Solution (HiCN) High-Normal level	2456	12 x 5 ml

SUMMARY

PRINCIPLES OF THE PROCEDURE

HiCN solution according the method described by Van Kampen en Zijlstra (1).

The hemiglobincyanide concentration, printed on the product label, is determined on a spectrophotometer based on a millimolar extinction coefficient of 44.0 for hemiglobincyanide at a wavelength of 540 nm and a molecule weight of 64,500 (2).

As reference preparation is used the Certified Reference Material (CRM) 552 and the 6th International Standard Haemiglobincyanide code 98/708 with a concentration of 49.8 micromoles/liter from the National institute for Biological Standards and Control, U.K.

PROCEDURE

Read the absorbance of the hemiglobincyanide solution on a spectrophotometer or filterphotometer at a wavelength of approx. 540 nm. Plot the absorbance (y-axis) versus hemoglobin concentrations (x-axis) on linear-linear graph paper. Use distilled water as a zero reference. The concentration of the measured hemoglobin concentration is known in $\mu\text{mol HiCN(Fe)/l}$. When a dilution factor of 251 is used, the corresponding hemoglobin concentration of the undiluted sample can be defined by multiply the result with 0.251.

OPTIONAL REAGENTS

Hemoglobin Reagent (Hemiglobincyanide)	2465	1 x 500 ml
Hemoglobin Reagent (Cyanide Free)	2488	1 x 500 ml

ANALYTICAL RESULTS FOR:

Article nr.: 2457, Hemiglobincyan. sol, low Contents: 12 x 5 ml

Component	Concentration	Component	Concentration
HiCN (Fe)	24.4 $\mu\text{mol/l}$	Hemoglobin	6.1 mmol/l

Article nr.: 2456, Hemiglobincyan. sol, High-normal Contents: 12 x 5 ml

Component	Concentration	Component	Concentration
HiCN (Fe)	35.3 $\mu\text{mol/l}$	Hemoglobin	8.9 mmol/l

NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstru[®]chemie for the complete product insert and latest edition.
4. Printed in the Netherlands, Hemiglobincyanide-Solutions-280910-1.FEN