



HEMOGLOBIN STANDARD

FOR USE AS A CALIBRATOR IN THE HEMOGLOBIN DETERMINATION

- Independent Calibrator
- Expiration date > 2 years, after production
- Ready for use
- Liquid



Products	Product no.	Quantity
Hemoglobin Standard Low Level (4.0 mmol/L)	2494	10 x 2 ml
Hemoglobin Standard Normal Level (8.2 mmol/L)	2495	10 x 2 ml
Hemoglobin Standard High Level (10.2 mmol/L)	2496	10 x 2 ml

SUMMARY

INTENDED USE

The Hemoglobin Standards (2494, 2495 and 2496) can be used to plot a hemoglobin calibration curve or as a single calibrator.

STABILITY

Hemoglobin Standards should be stored under refrigeration at 2-6 °C protected from extreme temperatures when in use.

OPTIONAL REAGENTS

Hemoglobin Reagent (Hemoglobin cyanide)	2465	1 x 500 ml
Hemoglobin Reagent (Cyanide Free)	2488	1 x 500 ml

NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstruChemie for the complete product insert and latest edition.
4. Printed in the Netherlands, Hemoglobin Standard -summary-280910-1.FEN