



## PHOSPHATE BUFFERED SALINE (INbusal)

- Concentrated PBS (10x)
- INbusal
- Expiration date > 8 years
- Store at room temperature



Products	Product no.	Quantity
Phosphate Buffered Saline - Concentrate (INbusal) PBS	2134	1 x 1 liter
Phosphate Buffered Saline - Container	2029	1 x 10 liter

## SUMMARY

### REAGENTS

INbusal Concentrate, product number 2134, is liquid. Store at room temperature.

Sodium Chloride

Disodium phosphate

Monosodium phosphate

For stabilities refer to expiry date printed on the product label.

### NOTES

1. For in vitro diagnostic use only.
2. For professional use only.
3. Always contact INstru<sup>®</sup>chemie for the complete product insert and latest edition.
4. Printed in the Netherlands, INbusal-PBS-summary-280901-1.FEN